

**IUE\***  
INTRA UTERINE ENABLERS



### **Intra Uterine Enablers (IUE)**

is a new series of breakthrough products researched and developed by Pregna International, a world leader in Intra-Uterine Contraceptive domain.

- \* IUE series of products focuses on uncompromising Client (women) safety, quick and easy procedure & enhanced ergonomics. IUE products are designed to follow Aseptic procedures, and reduce the anxious waiting time of Clients with an aim of making the entire experience for the providers and end users a safe and delightful one.
- \* State of art devices enabled to provide extra safety to end clients (Women)
- \* Intuitive devices enabling providers to use them with minimal training
- \* Feature enabled devices making the life of providers virtually carefree
- \* Technologically advanced products ensuring accuracy in use of products
- \* Designed with the end users as the starting point



*Introducing*

etherena 

An Integrated IUD Insertion Device



## Safe

Superior design offers perfect fundal placement thereby reducing bleeding, pain and possibility of expulsions

Reduces the possibility of accidental perforation of fundus

Works on principle of withdrawal technique only

Offers superior no touch technique, thus a very high aseptic level

Nylon suture is softer and hence is more comfortable for the couple

Device follows the natural curve of the vaginal/uterine cavity ensuring minimum discomfort for clients



## Intuitive

Completely integrated system for a safe and simple one hand operation

PP (Push-Pull) Simple two step process for quicker & easy insertion

Intuitive Click sound for knowing the set measurement for uterine depth



## Ergonomic

Device is operated by a single hand.

Ergonomic Design ensures that the service provider's fingers remain agile and fatigue free

Etherena comes with a special handle grip, making it easier for provider to handle and hold the device

Device comes with a curve similar to body's natural curve, thus allowing perfect insertion for women with any type of uterus (retroverted or anteverted)

Inbuilt uterine scale on the device for easy operability, convenience and accuracy



## Accurate

Uterine depth adjuster with precision up to 5mm

The loading of T-frame is accurate and requires minimum effort



Both the arms of T-frame are evenly loaded allowing correct insertion

---

Etherena has all the advantages of  
Copper T380A


# etherena

- \* No Touch Aseptic Device
- \* Enables accurate fundal placement of IUD
- \* Fatigue free operation
- \* Ergonomic handling
- \* Intuitive insertion process
- \* Quicker insertion/  
More patients attended





## 2 steps operation



Etherena makes the installation procedure of Copper-T 380A, a simple 2 step process. PP (Push-Pull) process makes the loading and inserting extremely intuitive.

A 2 step PP (Push-Pull) process deploys a withdrawal technique, ensuring an extremely safe and painless procedure, a first of its kind in the segment


# Know your product



Preloaded T held in place.  
No human contact required near the T during loading

Organic Curvilinear form  
adapts naturally to the  
uterine cavity ensuring  
accurate placement

Simple Push-Pull operation  
of the knob ensures quicker  
learning, faster error-free  
procedures





**13, Suryoday Estate, 136 Tardeo Road,  
Mumbai - 400034. India**

**Tel: +91 22 4345 4999**

**Fax: +91 22 2351 3484**

**Email: [sales@pregna.com](mailto:sales@pregna.com)**

**[www.pregna.com](http://www.pregna.com)**